

Determination of Enzyme Thermal Stability

The thermal stability of the enzyme was evaluated using a fluorescence dye-based Protein Thermal Shift (PTS) assay. This method monitors the changes in fluorescence intensity caused by the binding of a hydrophobic dye to the hydrophobic regions of proteins that become exposed during thermal denaturation. The experiments were conducted on an Applied Biosystems real-time PCR instrument (Thermo Fisher Scientific), with excitation and emission wavelengths set to 520 ± 10 nm and 558 ± 12 nm, respectively.

For each assay, 5 μ g of purified enzyme was mixed with 2.5 μ L of 8 \times Protein Thermal Shift dye. Subsequently, 5 μ L of reaction buffer was added, and the mixture was adjusted to a final volume of 20 μ L with phosphate buffer (50 mM Na₂HPO₄, 100 mM NaCl, pH 7.5). The prepared reaction mixture was transferred into the PCR instrument, equilibrated at 25 °C for 2 min, and then subjected to a temperature gradient from 25 °C to 99 °C with a ramp rate of 0.05 °C/min, while fluorescence was recorded in real time.

The fluorescence intensity profile reflected the thermal unfolding process of the protein. The melting temperature (T_m) was determined from the first derivative of the fluorescence versus temperature curve and used as a quantitative measure of enzyme thermal stability.